

Product Datasheet

Human HCC-1 protein (orb107735)

Description	Human CCL14 is belonging to the CC chemokine family. It is encoded by the gene CCL14. CCL14 has two ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	HCC-1
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 2xPBS, pH 7.4, 5 % trehalose
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	GPYHPSECCF TYTTYKIPRQ RIMDYETNS QCSKPGIVFI TKRGHSVCTN PSDKWVQDYIKDMKEN
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	7.8 kDa
Expiration Date	6 months from date of receipt.