


## Product Datasheet

### Porcine IL-1 beta protein (orb107741)

<b>Description</b>	Interleukin-1 beta (IL-1C) is a non-secreted proinflammatory cytokine produced mainly by activated m...
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	HPLC, SDS-PAGE
<b>Target</b>	IL-1 beta
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, 3 % trehalose
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	ANVQSMECKL QDKDHKSLVL AGPHMLKALH LLTGDLKREV VFCMSFVQGDDSNKIPVTL GIKGKNLYLS CVMKDNTPTL QLEDIDPKRY PKRDMEKRFVFKTEIKNRV EFESALYPNW YISTSQAEQK PVFLGNSKGR QDITDFTMEV LSP
<b>Purity</b>	> 95 % by SDS-PAGE and HPLC analyses
<b>MW</b>	17.6 kDa
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Porcine IL-1 beta p...