

## Product Datasheet

### IL-5 protein (orb107749)

**Description** Recombinant Rat Interleukin-5. IL-5, also named B-cell differentiation factor 1, eosinophil differen...

**Conjugation** Unconjugated


**Tested Applications** HPLC, SDS-PAGE

**Target** IL-5

**Storage** Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

**Note** For research use only

**Expiration Date** 6 months from date of receipt.



SDS-PAGE analysis of IL-5 protein