


## Product Datasheet

### Human IL-9 protein (orb107753)

<b>Description</b>	Interleukin-9 (IL-9) is encoded by the IL9 gene and produced by T-cells and specifically by CD4+ hel...
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	HPLC, SDS-PAGE
<b>Target</b>	IL-9
<b>Preservatives</b>	Lyophilized from a 0.2 m filtered concentrated solution in PBS, pH 7.4.
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	QGCPTLAGIL DINFLINKMQ EDPASKCHCS ANVTSCCLCG IPSDNCTRPC FSERLSQMTNTTMQTRYPLI FSRVKKSVEV LKNNKCPYFS CEQPCNQTTA GNALTFLKSL LEIFQKEKMRGMRGKI
<b>Purity</b>	> 95 % by SDS-PAGE and HPLC analyses
<b>MW</b>	14.1 kDa
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Human IL-9 protein