

Product Datasheet

Human IL-16(130a.a.) protein (orb107756)

Description	Human Interleukin-16 (IL-16) is also called lymphocyte chemoattractant factor (LCF) and it has three...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-16
Preservatives	Lyophilized from a 0.2 m filtered concentrated solution in 20 mM Tris, 150 mM NaCl, pH 7.0.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MPDLNSSTDS AASASAASDV SVESTAEATV CTVTLEKMSA GLGFSLEGGK GSLHGDKPLTINRIFKGAAS EQSETVQPGD EILQLGGTAM QGLTRFEAWN IIKALPDGPV TIVIRRKSLQSKETTAAGDS
Purity	> 95 % by SDS-PAGE and HPLC analyses.
MW	13.4 kDa
Expiration Date	6 months from date of receipt.