


Product Datasheet

Murine IL-22 protein (orb107759)

Description	IL-22 is belonging to IL-10 family of regulatory cytokines which includes IL-10, IL-19, IL-20, IL-22...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-22
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	LPVNTRCKLE VSNFQQPYIV NRTFMLAKEA SLADNNTDVR LIGEKLFGRGVSADQCYLMK QVLNFTLEDV LLPQSDRFQP YMQEVPFLT KLSNQLSSCHISGDDQNIQK NVRRLKETVK KLGESGEIKA IGELDLLFMS LRNACV
Purity	> 96 % by SDS-PAGE and HPLC analyses.
MW	16.6 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Murine IL-22 protei...