

Product Datasheet

Human IL-36 alpha (153a.a.) protein (orb107763)

Description	Interleukin-36 (IL-36) is a pro-inflammatory cytokine which plays an important role in the pathophys...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-36 alpha
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris, 300 mM NaCl, pH 8.0, 0.1 % Tween 80.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	KIDTPQQGSI QDINHRVWVL QDQTLIAVPR KDRMSPVTIA LISCRHVETL EKDRGNPIYLGLNGLNLCLM CAKVGDQPTL QLKEKDIMDL YNQPEPVKSF LFYHSQSGRN STFESVAFPGWFIASSEGG CPLILTQELG KANTTDFGLT MLF
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	17.1 kDa
Expiration Date	6 months from date of receipt.