


## Product Datasheet

### Murine IL-36 alpha (160a.a.) protein (orb107764)

<b>Description</b>	Interleukin-36 (IL-36) is a pro-inflammatory cytokine which plays an important role in the pathophys...
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	HPLC, SDS-PAGE
<b>Target</b>	IL-36 alpha
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4, 5 % trehalose.
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MNKEKELRAA SPSLRHVQDL SSRVWILQNN ILTAVPRKEQ TVPVTITLLP CQYLDTLETNRGDPTYMGVQ RPMSCFCTK DGEQPVLQLG EGNIMEMYNK KEPVKASLFYHKKSGTTSTF ESAAFPGWFI AVCSKGSCPL ILTQELGEIF ITDFEMIVVH
<b>Purity</b>	> 95 % by SDS-PAGE and HPLC analyses
<b>MW</b>	18.0 kDa
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Murine IL-36 alpha ...