


Product Datasheet

Murine IL-36 beta (183a.a.) protein (orb107767)

Description	Interleukin-36 (IL-36) is a pro-inflammatory cytokine which plays an important role in the pathophys...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-36 beta
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris, 300 mM NaCl, pH 8.0, 5 %trehalose.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MMAFPPQSCV HVLPPKSIQM WEPNHNTMHG SSQSPRNYRV HDSQQMVWVLTGNTLTAVPA SNNVKPVILS LIACRDTEFQ DVKKGNLVFL GIKNRNLCFCCVEMEGKPTL QLKEVDIMNL YKERKAQKAF LFYHGIEGST SVFQSVLYPGWFIATSSIER QTIILTHQRG KLVNTNFYIE SEK
Purity	> 95 % by SDS-PAGE and HPLC analyses
MW	20.9 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Murine IL-36 beta (...)