

Product Datasheet

Human IL-36 gamma (169a.a.) protein (orb107768)

Description	Interleukin-36 (IL-36) is a pro-inflammatory cytokine which plays an important role in the pathophys...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	IL-36 gamma
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MRGTPGDADG GGRAVYQSMC KPITGTINDL NQQVWTLQGQ NLVAVPRSDSVTPVTAVIT CKYPEALEQG RGDPIYLGIQ NPEMCLYCEK VGEQPTLQLK EQKIMDLYGQPEPVKPFIFY RAKTGRSTL ESVAFPDWFI ASSKRDQPII LTSELGKSYN TAFELNIND
Purity	> 95 % by SDS-PAGE and HPLC analyses.
MW	18.7 kDa
Expiration Date	6 months from date of receipt.