

## Product Datasheet

### Murine IP10 protein (orb107772)

<b>Description</b>	Recombinant Murine gamma-Interferon Inducible Protein 10/CXCL10. CXCL10 also known as IP-10 is belon...
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	HPLC, SDS-PAGE
<b>Target</b>	IP10
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	IPLARTVRCN CIHIDDPVVR MRAIGKLEII PASLSCPVE IIATMKKNDK QRCLNPESKTIKMLKAFSQ KRSKRAP
<b>Purity</b>	> 97% by SDS-PAGE and HPLC analyses.
<b>MW</b>	8.7 kDa
<b>Expiration Date</b>	6 months from date of receipt.