

Product Datasheet

Human Neuroserpin protein (orb107789)

Description	Neuroserpin is an inhibitory serpin that is expressed predominantly in central nervous system. Altho...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Neuroserpin
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MTGATFPPEA IADLSVNMYN RLRATGEDEN ILFSPLSIAL AMGMMELGAQGSTQKEIRHS MGYDSLKNGE EFSFLKEFSN MVTAKESQYV MKIANSLFVQNGFHVNEEFL QMMKKYFNAA VNHVDFSQNV AVANYINKKW ENNTNNLVKDLVSPRDFDAA TYLALINAVY FKGNWKSQFR PENTRTFSFT KDDESEVQIPMMYQQGEFY GEFSDGSNEA GGIYQVLEIP YEGDEISMML VLSRQEVPLATLEPLVKAQL VEEWANSVKK QKVEVYLPRF TVEQEIDLKD VLKALGITEIFIKDANLTGL SDNKEIFLSK AIHKSFLEVN EEGSEAAAVS GMIAISRMVLYPQVIVDHP FFFLIRNRRT GTILFMGRVM HPETMNTSGH DFEEL
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	44.8kDa
Expiration Date	6 months from date of receipt.