

Product Datasheet

Murine PF-4 protein (orb107796)

Description	Platelet Factor 4, is a member of the CXC chemokine family, CXCL4. CXCL4 has homology with IL8 and b...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	PF-4
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in 20mM PB, 1.5M NaCl, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	VTSAGPEESD GDLSCVCVKT ISSGIHLKHI TSLEVIKAGR HCAVPQLIAT LKNRKRICLDRQAPLYKKVI KKILES
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	8.2 kDa
Expiration Date	6 months from date of receipt.