

Product Datasheet

Human PTH(1-84) protein (orb107797)

Description	Human PTH(1-84) protein
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Human PTH(1-84)
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	SVSEIQLMHN LGKHLNSMER VEWLRKKLQD VHNFVALGAP LAPRDAGSQRPRKKEDNVLV ESHEKSLGEA DKADVNVLTAKSQ
Purity	> 97 % by SDS-PAGE and HPLC analyses
MW	9.5 kDa
Expiration Date	6 months from date of receipt.