


Product Datasheet

Murine SCF protein (orb107800)

Description	C-kit ligand, the recently identified ligand for the kit tyrosine kinase receptor, is mapped to the ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	SCF
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MKEICGNPVT DNVKDITKLV ANLPNDYMIT LNYVAGMDVL PSHCWLRDMVIQLSLSLTTL LDKFSNISEG LSNYSIIDKL GKIVDDL VLC MEENAPKNIK ESPKRPETRSFTPEEFSIF NRSIDAFKDF MVASDTSDCV LSSTLGPEKD SRVSVTKPFM LPPVA
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	18.4 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Murine SCF protein