


Product Datasheet

Streptokinase protein (orb107804)

Description	Streptokinase is a thrombolytic agent i.e. it dissolves vascular thrombi. This potent agent is deriv...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	Streptokinase
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	IAGPEWLLDR PSVNNSQLVV SVAGTVEGTN QDISLKFFEI DLTSRPAHGG KTEQGLSPKS KLFATDSGAM PHKLEKADLL KAIQEQLIAN VHSNDDYFEV IDFASDATIT DRNGKVYF ADKDGSVTLPIQ PVQEFLKKGH VRVRPYKEKP VQNQAKSVDV EYTVQFTPLN PDDDFRPA LKDTKLLKTLAI GDTITSQELL AQAQSILNKN HPGYTIYERD SSIVTHDNDI FRTILPMDQE FTYHVKNREQ AYRINKKSGL NEEINNTDLI SEKYYVLKKG EKPYPDFDRS HCLKFTIKYV DVNTNELLKS EQLLTASERN LDFRDLYDPR DKAKLLYNNL DAFGIMDYT LTGKVEDNHDD TNRIITVYMG KRPEGENASY HLAYDKDRYT EEEREVYSY LRYTGTPIDN PNDK
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	47.3 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Streptokinase prote...