


## Product Datasheet

### Murine TARC protein (orb107806)

|                            |                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Thymus and activation regulated chemokine (TARC) is a novel CC chemokine, also called CCL17, recentl...        |
| <b>Conjugation</b>         | Unconjugated                                                                                                   |
| <b>Tested Applications</b> | HPLC, SDS-PAGE                                                                                                 |
| <b>Target</b>              | TARC                                                                                                           |
| <b>Preservatives</b>       | Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.                                       |
| <b>Storage</b>             | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                                          |
| <b>Protein Sequence</b>    | ARATNVGREC CLDYFKGAIP IRKLVSWYKT SVECSRDAIV<br>FLTVQGKLIC ADPKDKHVKKAIRLVKNRP                                  |
| <b>Purity</b>              | > 97% by SDS-PAGE and HPLC analyses                                                                            |
| <b>MW</b>                  | 7.9kDa                                                                                                         |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                 |



SDS-PAGE  
analysis of  
Murine TARC  
protein