


Product Datasheet

Human UBC2 protein (orb107812)

Description	This E2 enzyme encodes for the human homolog of the yeast DNA repair gene RAD6, which is induced by ...
Conjugation	Unconjugated
Tested Applications	HPLC, SDS-PAGE
Target	UBC2
Preservatives	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MHHHHHHAMG QLRSMSTPAR RRLMRDFKRL QEDPPVGVSG APSENNIMQWNAVIFGPEGT PFEDGTFKLV IEFSEEYPNK PPTVRFLSKM FHPNVYADGS ICLDILQNRWSPTYDVSSIL TSIQSLLDEP NPNSPANSQA AQLYQENKRE YEKRVSAIVE QSWNDS
Purity	> 95% by SDS-PAGE and HPLC analyses
MW	19.0 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human UBC2 protein