


Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody (orb1084018)

Description	SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 (2019-nCoV) Nucleocapsid protein (full length)
Target	SARS-CoV-2 (2019-nCoV) Nucleocapsid
Preservatives	0.01M PBS(pH7.4) with 0.03% Proclin300
Form/Appearance	Liquid
Concentration	≥1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	PODTC2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis using Anti-SARS-Co...