

Product Datasheet

ADGRE1 (F4/80) antibody (orb1084093)

| | |
|----------------------------|---|
| Description | ADGRE1 (F4/80) antibody |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | KLH conjugated synthetic peptide derived from mouse F4/80 |
| Target | ADGRE1 (F4/80) |
| Preservatives | 0.01M PBS(pH7.4) with 0.1% rAlbumin, 0.03% Proclin300. |
| Form/Appearance | Liquid |
| Concentration | 200ug/ml |
| Storage | Shipped at 4°C. Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Monoclonal |
| MW | 95 kDa |
| Uniprot ID | Q14246 |
| Dilution Range | WB=1:100-1000, IHC-P=1:100-500, IF=1:50-500 |
| Expiration Date | 12 months from date of receipt. |