

Product Datasheet

cyclin E antibody (orb108539)

Description	Mouse monoclonal to cyclin E
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant full-length human cyclin E protein
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody detects human cyclin E protein by Western blot. Detects a 50 kDa band. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P24864
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.