

Product Datasheet

Human PBEF/Visfatin ELISA KIT (orb1085892)

Description	Human PBEF/Visfatin ELISA KIT
Reactivity	Human
Range	3.125 - 200 ng/mL
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum, and plasma.
Sensitivity	2 ng/mL
Expiration Date	Please enquire.