

## Product Datasheet

### Human lipocalin-2/NGAL ELISA Kit (orb1085930)

<b>Description</b>	Human lipocalin-2/NGAL ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78 - 5,000pg/mL
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).
<b>Note</b>	For research use only
<b>Sample Types</b>	Cell culture supernates, serum, and plasma.
<b>Sensitivity</b>	40 pg/mL
<b>Expiration Date</b>	Please enquire.

---

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)