


Product Datasheet

p107 antibody (orb108603)

Description	KAB6 recognizes both the phosphorylated and unphosphorylated forms of p107 at a MW of 107kDa.
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Western Blots at 1:1000
Isotype	IgG1
Clonality	Monoclonal
Source	The immunogen used to immunize the mice was a recombinant GST protein coupled to the N-terminus of p107.
Uniprot ID	Q64701, P28749
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant GST...