

Product Datasheet

Human HSP70/HSPA1A ELISA KIT (orb1086132)

Description	Human HSP70/HSPA1A ELISA KIT
Reactivity	Human
Range	125.00 - 8,000 pg/mL
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum, and plasma.
Sensitivity	60 pg/mL
Expiration Date	Please enquire.