

Product Datasheet

Human CXCL11/I-TAC ELISA KIT (orb1086144)

Description	Human CXCL11/I-TAC ELISA KIT
Reactivity	Human
Range	7.81 - 500 pg/mL
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum, and plasma.
Sensitivity	4 pg/mL
Expiration Date	Please enquire.