

## Product Datasheet

### Human CXCL6/GCP-2 ELISA KIT (orb1086145)

<b>Description</b>	Human CXCL6/GCP-2 ELISA KIT
<b>Reactivity</b>	Human
<b>Range</b>	15.60 - 1,000 pg/mL
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernates, serum, and plasma.
<b>Sensitivity</b>	7 pg/mL
<b>Expiration Date</b>	Please enquire.