

## Product Datasheet

### Human CCL26/Eotaxin-3 ELISA KIT (orb1086149)

<b>Description</b>	Human CCL26/Eotaxin-3 ELISA KIT
<b>Reactivity</b>	Human
<b>Range</b>	62.50 - 4,000 pg/mL
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernates, serum, and plasma.
<b>Sensitivity</b>	31 pg/mL
<b>Expiration Date</b>	Please enquire.