

Product Datasheet

Human CXCL2/GRO β ELISA KIT (orb1086186)

Description	Human CXCL2/GRO β ELISA KIT
Reactivity	Human
Range	15.625 - 1000 pg/mL
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum, and plasma.
Sensitivity	7 pg/mL
Expiration Date	Please enquire.