

Product Datasheet

Mouse IL-25/IL-17E ELISA KIT (orb1086231)

Description	Mouse IL-25/IL-17E ELISA KIT
Reactivity	Mouse
Range	31.25-2000 pg/mL
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum, and plasma.
Sensitivity	16 pg/mL
Expiration Date	Please enquire.