

Product Datasheet

Mouse CXCL9/MIG ELISA KIT (orb1086264)

| | |
|------------------------|--|
| Description | Mouse CXCL9/MIG ELISA KIT |
| Reactivity | Mouse |
| Range | 15.60-1000 pg/mL |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.) |
| Note | For research use only |
| Sample Types | cell culture supernates, serum, and plasma. |
| Sensitivity | 7 pg/mL |
| Expiration Date | Please enquire. |