

Product Datasheet

Mouse BAFF/BLyS/TNFSF13B ELISA KIT (orb1086269)

Description	Mouse BAFF/BLyS/TNFSF13B ELISA KIT
Reactivity	Mouse
Range	62.50-4000 pg/mL
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum, and plasma.
Sensitivity	30 pg/mL
Expiration Date	Please enquire.