

Product Datasheet

Mouse CCL24/Eotaxin/MPIF-2 ELISA KIT (orb1086278)

Description	Mouse CCL24/Eotaxin/MPIF-2 ELISA KIT
Reactivity	Mouse
Range	15.60 - 1,000 pg/mL
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.)
Note	For research use only
Sample Types	cell culture supernates, serum, and plasma.
Sensitivity	7 pg/mL
Expiration Date	Please enquire.