

Product Datasheet

Protein Phosphatase V/C antibody (orb108639)

Description	Sheep polyclonal to Protein Phosphatase V/C
Species/Host	Sheep
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A peptide conjugated to KLH
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Applications: Immunoblotting 1:500 to 1:1,000. Immunoprecipitation - Empirical.
Isotype	IgG
Uniprot ID	Q9CQR6 , O00743
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.