

Product Datasheet

Mouse AGR2 / AG-2 Protein (orb1087549)

Description	Mouse AGR2, His Tag is expressed from human 293 cells (HEK293). It contains AA Lys 21 - Leu 175 (Acc...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	AGR2
Preservatives	PBS, 1 mM TCEP, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 19.8 kDa. The protein migrates as 23 kDa and 24 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_035913.1
Purity	95%
Source	Mouse AGR2, His Tag (orb1087549) is expressed from human 293 cells (HEK293). It contains AA Lys 21 - Leu 175 (Accession # O88312-1).
MW	19.8 kDa
Uniprot ID	O88312
NCBI	NP_035913.1
Expiration Date	6 months from date of receipt.