

## Product Datasheet

### Human GAS6 (279-678) Protein (orb1087554)

<b>Description</b>	Human GAS6, His Tag is expressed from human 293 cells (HEK293). It contains AA Asp 279 - Ala 678 (Ac...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	GAS6
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	N-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 46.0 kDa. The protein migrates as 45-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NP_000811.1
<b>Purity</b>	95%
<b>Source</b>	Human GAS6, His Tag (orb1087554) is expressed from human 293 cells (HEK293). It contains AA Asp 279 - Ala 678 (Accession # Q14393-2).
<b>MW</b>	46.0 kDa
<b>Uniprot ID</b>	<a href="#">Q14393</a>
<b>NCBI</b>	<a href="#">NP_000811.1</a>
<b>Expiration Date</b>	6 months from date of receipt.