

Product Datasheet

Enzyme (orb1087603)

Description	Enzyme
Tested Applications	IF, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Storage	The solution should be stored at 2-8°C without further dilution. Under these conditions the solution is stable up to the expiry date indicated on the label. Do not use product after the expiry date. A positive and a negative control have to be carried out in parallel to the test material. If you observe unusual staining or other deviations from the expected results which could possibly be caused by this reagent, please contact our technical support.
Note	For research use only
Application notes	Fast Enzyme is a ready-to-use enzyme solution for enzymatic epitope retrieval.
Expiration Date	12 months from date of receipt.