

## Product Datasheet

# Mouse anti-IgG2a Isotype control, clone MOPC-173, Biotin (Monoclonal) (orb1088085)

<b>Description</b>	Mouse anti-IgG2a Isotype control, clone MOPC-173, Biotin (Monoclonal)
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC, ICC, IHC-Fr, IHC-P, IP, WB
<b>Immunogen</b>	Please contact us for more information.
<b>Preservatives</b>	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a kappa
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	.....
<b>Expiration Date</b>	12 months from date of receipt.