

Product Datasheet

Sheep SELP(Selectin, Platelet) ELISA Kit (orb1088126)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Sheep
Range	6.25-400 ng/mL
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P98109
Sensitivity	2.73 ng/mL
Expiration Date	Please enquire.

