

## Product Datasheet

# Human RAD51B(DNA repair protein RAD51 homolog 2) ELISA Kit (orb1088141)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 0.32-20 ng/mL

**Note** For research use only

**Sample Types** Tissue homogenates, cell lysates and other biological fluids.

**Uniprot ID** [O15315](#)

**Sensitivity** 0.115 ng/mL

**Expiration Date** Please enquire.

