

Product Datasheet

Human PIK3C2a(Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit (orb1088144)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Note	For research use only
Sample Types	Tissue homogenates and other biological fluids.
Uniprot ID	O00443
Sensitivity	0.057 ng/mL
Expiration Date	Please enquire.

