

Product Datasheet

Human sCTLA4(Soluble Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit (orb1088146)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	15.63-1000 pg/mL
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids
Sensitivity	9.38 pg/mL
Expiration Date	Please enquire.

