

## Product Datasheet

### Human Hepc(Hepcidin) ELISA Kit (orb1088150)

|                        |   |
|------------------------|---|
| <b>Description</b>     | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| <b>Reactivity</b>      | Human   |
| <b>Range</b>           | 0.63-40 ng/mL   |
| <b>Note</b>            | For research use only   |
| <b>Sample Types</b>    | Tissue homogenates and other biological fluids.   |
| <b>Uniprot ID</b>      | <b>P81172</b>   |
| <b>Sensitivity</b>     | 0.218 ng/mL   |
| <b>Expiration Date</b> | Please enquire.   |

