

## Product Datasheet

# Human PGH2(Prostaglandin H2) ELISA Kit (orb1088173)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Human
<b>Range</b>	3.13-200 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Sensitivity</b>	0.97 pg/mL
<b>Expiration Date</b>	Please enquire.

