

## Product Datasheet

### Zebrafish DA(Dopamine) ELISA Kit (orb1088183)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Zebrafish
<b>Range</b>	15.63-1000 pg/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma and other biological fluids
<b>Sensitivity</b>	4.71 pg/mL
<b>Expiration Date</b>	Please enquire.

