

## Product Datasheet

# Rat ATP(Adenosine Triphosphate) ELISA Kit (orb1088186)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Rat
<b>Range</b>	15.63-1000 ng/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Sensitivity</b>	5.27 ng/mL
<b>Expiration Date</b>	Please enquire.

