

Product Datasheet

Human Zinc transporter ZIP8 (SLC39A8) ELISA Kit (orb1088213)

| | |
|------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 0.16-10 ng/mL |
| Note | For research use only |
| Sample Types | tissue homogenates, cell lysates and other biological fluids |
| Uniprot ID | Q9C0K1 |
| Sensitivity | 0.065 ng/mL |
| Expiration Date | Please enquire. |

