


## Product Datasheet

# Human PC(Protein Carbonyl) ELISA Kit (orb1088214)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Human
<b>Range</b>	0.32-20 nmol/mL
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Cell Lysate, Hep Plasma, EDTA Plasma, Tissue Lysate
<b>Sensitivity</b>	0.119 nmol/mL
<b>Expiration Date</b>	Please enquire.



Standard Curve of Human PC(Protein Carbonyl) ELISA Kit

