
Product Datasheet

CD172g Antibody (PE) (orb1088476)

Description

CD172g is a transmembrane glycoprotein, which may

Reactivity

Human

Conjugation

PE

Tested Applications

FC

Immunogen

recombinant human CD172g

Target

CD172g

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

[Q9P1W8](#)

Entrez


[55423](#)

Dilution Range


Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Expiration Date

12 months from date of receipt.



Separation of human CD172g positive lymph...



Flow cytometry surface staining pattern ...