



Product Datasheet

HLA-C Antibody (orb1088498)

Description	HLA-C, a member of MHC class I glycoproteins, is one of polymorphysm typing targets, which are impor...
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Purified MHC class I molecules of tamarin origin
Target	HLA-C
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P10321
Entrez	3107
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of lymphocytes stained anti-h...



Flow cytometry surface staining pattern ...